

PSK ENDO-SUCTION PRE-CLEANING SYSTEM

PSK USAGE TRACKING SHEET*

***Maximum of 500 uses per PSK unit**



PSK Model No. _____

PSK Product Serial No. _____

Tel: +61 3 9887 7629 www.pskendoscopy.com.au

Fax: +61 3 9886 1888 info@pskendoscopy.com.au

Instructions

1. Fill in the table below at the start of each '25 PSK usage' cycle
2. Start using the PSK unit to pre-clean your scope
3. After you reach 25 uses, photocopy this tracking sheet and re-start the process until you reach 450 PSK uses.
4. Re-order a new replacement PSK unit so it arrives before you reach 500 uses

Start Date of 25 use cycle: ____ / ____ / ____

| Component/Item to Check | Colour Change | | Cracks | |
|--|--------------------------|--------------------------|--------------------------|--------------------------|
| | No change - clear tubing | Change | None | Requires Attention |
| BLUE Tubing marked Air/Water | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| RED Tubing marked Suction | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| RED Tubing marked Biopsy | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| White Tubing marked PSK | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| PSK Suction Cap (For Olympus Endoscopes) | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| PSK Air/Water Cap (For Olympus Endoscopes) | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |

TOTAL USAGE FOR THIS SERIAL NO. _____ / 500 maximum

More copies available on www.pskendoscopy.com.au/service-support/page-2.html

You may also make copies of this tracking sheet.

| Daily Procedure No. | Endoscope Serial No. | Date | Sign-off Signature |
|---------------------|----------------------|------|--------------------|
| 1 | | / / | |
| 2 | | / / | |
| 3 | | / / | |
| 4 | | / / | |
| 5 | | / / | |
| 6 | | / / | |
| 7 | | / / | |
| 8 | | / / | |
| 9 | | / / | |
| 10 | | / / | |
| 11 | | / / | |
| 12 | | / / | |
| 13 | | / / | |
| 14 | | / / | |
| 15 | | / / | |
| 16 | | / / | |
| 17 | | / / | |
| 18 | | / / | |
| 19 | | / / | |
| 20 | | / / | |
| 21 | | / / | |
| 22 | | / / | |
| 23 | | / / | |
| 24 | | / / | |
| 25 | | / / | |