

SIMPLIFYING ENDOSCOPE PRE-CLEANING...

Introducing the PSK Endo-Suction Cleaning System (ESCS), a revolutionary and innovative way of quickly pre-cleaning your scopes. Using suction rather than the traditional “push” method, it has been used world-wide for over 16 years.

Ask for a demonstration today!

✓ MOST EFFICIENT AND EFFECTIVE PRE-CLEANING SYSTEM IN INDUSTRY

✓ MAY BE USED TO MANUALLY CLEAN ENDOSCOPES OR PRE-CLEAN ONLY

✓ THE COST OF AN ESCS IS LESS THAN A GASTROSCOPY PROCEDURE

✓ ONLY WORLD-WIDE SYSTEM TO USE SUCTION RATHER THAN “PUSH”

✓ FAST AND SIMPLE TO USE, WITH ALL MAJOR SCOPE BRANDS CATERED FOR


✓ REDUCES RSI BY BEING A SEMI-AUTOMATED SYSTEM

✓ AUTOCLAVABLE


✓ GUARANTEED UP TO 500 USES PER ESCS*

✓ CLEAR EVIDENCE OF VISUAL CLEANING (UNLIKE AUTOMATED PUSH SYSTEMS)


✓ NO COSTLY REGULAR MAINTENANCE OR CALIBRATION (UNLIKE AUTOMATED PUSH SYSTEMS)



Example of an ESCS for a Pentax Scope



Example of an ESCS for an Olympus Scope



Suction accessory equipment for the ESCS

* DEPENDENT ON USAGE AND DETERGENTS APPLIED

Manufactured by:
PSK Endoscopy Australia Pty Ltd
55 Marianne Way Mt. Waverley
VIC 3149 AUSTRALIA
☎ (03) 9887 7629
✉ (03) 9886 1888
✉ info@pskendoscopy.com.au
🌐 www.pskendoscopy.com.au



“The Experts in Cleaning Your
Scopes Smarter”

Proudly Distributed in WA by:
Western Biomedical
41 Oxford Close, West Leederville
WA 6007 AUSTRALIA
☎ (08) 9381 9944
✉ (08) 9388 1715
✉ admin@westernbiomedical.com.au
🌐 www.westernbiomedical.com.au

